

# DETERGENTS & CLEANER NU®



## NU® GLASS WASHING TABLETS

Perfekt hygiene!  
500g can with 100 glass detergent tablets for 100 applications.  
Guaranteed hygienically pure, beer foam-friendly glasses. Already suitably dosed, no quarter of the tablet necessary, optimum and dry storage in the can.



art-no. 1700410

## NU® BRUSH SANITIZER

Desinfection and care!  
750g can of brush sanitizer for about 25 applications. SPÜLBOY® brush sanitizer ensures constantly fresh brushes and 100% hygiene.



art-no. 1700412

## Easy to use!

Add 1 NU® glass detergent tablet into the pre-rinsing canister (brush canister) and, depending on the quantity of glass washes, rinse 100% hygienically for 6-8 hours.

**Hint!** If no foam forms, the tablet has completely dissolved. Time to refill another 1 tablet.

## Clear advantages!

- Particularly slow release of the NU® tablets guarantees 6-8 hours of **hygienic cleaning power**.
- Optimum activation and distribution of hygiene agents through the unique **circulation of fresh water**
- Gentle ingredients of the NU® tablets – **longer life** of the brushes
- Low foaming and **only!** during the active rinsing process. So you have a clear view for the handling of glasses.

## Better than liquid!

- **Optimal effect!** NU® tablets do not float on the water, making them last longer.
- **Sustainable hygiene!** NU® glass washing tablets are biodegradable and extremely economical.

## Daily Hygiene Plus! NU® Brush cleaner

To care for and disinfect each SPÜLBOY® glass cleaning brush: Carefully remove the brushes, cover with 40°C warm water, add **2 parts** of clean brush and leave for around 8 hours (overnight). Rinse well in the morning and you're brushes are clean an in a good shape again!